Italiano

Spagnolo

Inglese

Italiano

Biofilie medium

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English

Biofilie medium

We created this recipe with the intention of having a medium for bacterial cultivation that is easy to make at home. We tested this protocol in Barcelona at Hangar Wetlab, retrieving all the ingredients without buying them on the Internet and going to the supermarket. Therefore, keep in mind that their cost and availability is highly subjective to where we worked and we will be interested in understanding the accessibility of ingredients in different locations. We then cultivated both bacteria already isolated and those present in our vaginal fluid, with excellent results. Below you will find all the information you need to carry out the preparation, but if you have any doubts, do not hesitate to ask us. We would then be curious to receive your feedback and if necessary change the protocol.

The Biofilie Lab project is part of a larger work called the Alma Toolkit, a set of educational tools and practices for self-exploration, a fundamental means carried out by feminist collectives to reappropriate the body. Microbiology, science and laboratories are still inaccessible but fundamental fields because they are not neutral and in turn create categories and identities that are subject to a hierarchy of knowledge. This is why we decided to develop this protocol and focus on the vaginal microbiome, i.e. the set of micro-organisms with which we live in symbiosis. We hope that this protocol will become more and more accessible and that it will serve to better understand this symbiosis.

About the vaginal microbiota

The microbiota is the collection of micro-organisms present in, on and with the human body. It includes bacteria, but also viruses, mycetes and protozoa. Science has subdivided it into districts located in different areas of the body, so we speak of vaginal, intestinal, oral, skin, etc. microbiota, each one different from the other but actually in close communication and symbiosis. In the vaginal area, there is a homeostatic and mutualistic relationship between the microbiota and the human host: the host provides a moist, nutritious and warm habitat for the microbes, while the resident microbiota produces antimicrobial and anti-inflammatory factors. Most vaginal microbiomes are dominated by single or multiple *Lactobacillus* species from *L. crispatus*, *L. gasseri*, *L. iners* and *L. jensenii* respectively, while some people may have a high diversity of microbial community characterised by other anaerobic bacteria.

Lactobacillus species thrive in the anaerobic vaginal environment and produce lactic acid to maintain a low pH level (between 3.8 and 4.5), which helps protect the vaginal microbiome from harmful pathogens. They occupy space on the vaginal wall to prevent other harmful bacteria from remaining there and produce bacteriocins (the body's natural antibiotic) to prevent the growth of harmful bacteria. When there are not enough protective bacteria such as *Lactobacilli* in the vaginal microbiome, disturbing bacteria such as *Gardnerella*, *E.coli* and *Mycoplasma* can take over and cause vaginal dysbiosis. Vaginal dysbiosis is a generic term used to describe a state of imbalance of the vaginal microbiome. The most common type of

vaginal dysbiosis is an infection called bacterial vaginosis (BV), which is clinically defined as an overgrowth of disruptive bacteria in the vaginal microbiome.

The dominant *Lactobacillus* species determines the extent to which the vaginal ecosystem is protected. For instance, dysbiosis and low stability are usually correlated with a vaginal microbiota dominated by *L. iners*. In contrast, the health and high stability of the vaginal community is enhanced by *L. crispatus*, which provides the necessary lactic acids.

About bacterial cultivation

Bacterial cultivation is a laboratory technique used to identify bacterial colonies in our vaginal fluid. In fact, a specific medium is created in which there are foods that the bacteria are fond of, thus promoting their growth so that they can be observed. Nowadays techniques are becoming more and more sophisticated, but we have decided to recover this kind of knowledge and practice because we consider it accessible as a tool to get to know our bodies.

Although it is also used for the diagnosis of specific infections, in our case the protocol is mostly designed to get in touch with the vaginal microbiome and open up reflections and discussions.

The Biofilie protocol works well with fungi (*Candida albicans*) and lactobacilli (*Gasseri* and *Jensenii*). It would also work with *Gardnerella*, but being an anaerobic bacterium, which therefore does not survive in contact with oxygen, we found it difficult to develop the technique for this bacterium.

Ingredients and costs (in euros and approximate):

- Fish food. Price 4 euro
- Polysorbate 80: This is a hyposoluble fat used to make soaps. We found it available in a shop but it is much more likely that you have to buy it on the internet. If you can think of other water-soluble fats, suggestions are welcome! Price: 6 euro
- Yeast extract. Commercial name is marmite but as we could not find it we used Viandel. Price: 5 euro (marmite on the internet costs 15 euro)
- AgarAgar (vegetable gelatine). Price 8 euro
- Distilled water. Price: 3 euro

Total cost: 25-30euro

Equipment

- Pressure cooker: an ordinary pressure cooker that works well in place of an autoclave for sterilising specific media. Price: 30 euros
- Borosilicate glass bottle: a borosilicate glass bottle where you add the ingredients and insert into the pressure cooker to sterilise the contents. Price: 5-10euro
- Autoclave tape: an adhesive tape that is coloured to confirm that sterilisation has taken place. It is not obligatory but we recommend it if accessible. Price: 10euro

- Inoculators: we used metal blackhead removers :) which work really well as inoculators and are easy to sterilise and reuse. Price: 5-9euros
- Plastic or glass Petri dishes: essential for bacterial culture. Price: 0.70 cent per petri (plastic)
- Parafilm: is a type of elastic, waxed tape with which you can seal the petri dish to prevent contamination but allow steam to escape. Price: 3 euro
- Candles: necessary for when you pour the medium into the petri dish and for inoculating the bacteria if you do not have a sterilisation chamber. Price 10 euro (x100 pieces)
- Swabs: for collecting vaginal fluid. They look like cotton buds but longer. Price: 9 euros
- Speculum: can be found in pharmacies, although from experience we have found that in some places you have to insist a bit to get it. Price: 1 euro
- Vaginal tampon: can be found in pharmacies. Price: 1.10 euro
- Egg incubator: works very well for incubating petri dishes and allowing bacteria to grow. You can use either a chicken egg incubator or a snake egg incubator. Price: 30 euro
- Mirror and light to observe your cervix: we recommend a small, portable mirror and a portable book light or phone flash

Total cost of equipment: 133 euro

Instructions

*Contamination is an interesting phenomenon but it is best to avoid it to see if the soil is working and candida or lactobacilli are growing. So always take care that the instruments you use are sterile, use 96% alcohol diluted with water, gloves, etc..

Prepare the medium (considering 1l of water)

15 g powdered agar

10gr fish food

*Tip: pulverise it and infuse it in distilled water at 60°. Remember that distilled water takes a short time to evaporate. Filter and add to the mixture.

4 g yeast extract

1gr Polysorbate 80

2gr of 4-5% white vinegar

1l distilled water

Mix well!

Place in the pressure cooker, with 20g agave on the side to prevent the mixture from caramelising. From when the pressure cooker starts to whistle count 20minutes. To make sure it has sterilised, we recommend using the scotch tape provided for this purpose, which will turn from green to black if sterilisation has taken place.

Mix the agave with the rest of the mixture.

Preparing petri dishes

In a laboratory, a sterilisation chamber is usually used for both petri dish preparation and inoculation. Or you can build a sterilisation chest but we were not happy with this and preferred to use the ancient method of candles. Place five candles as if to form a five-pointed star, work within that area, trying not to talk and breathe through your nose, sprinkle alcohol every now and then.

Before the medium hardens, pour the contents into the petri dishes inside the star of candles.

As you prepare them, seal the petri dishes with parafilm and place them in a box with candles inside to allow them to cool until they harden (time about 1 hour)

Collecting the sample

You can collect the sample without needing the speculum, but from experience we have seen that more fluid can be collected this way with less epithelial cells.

Take the speculum and use it as shown in the picture

Once the cervix is found, insert the swab and collect fluid between the cervix and the vaginal canal

Lactobacilli are semi-anaerobic in the sense that they have some tolerance to oxygen but do not survive long so do not wait too long to inoculate them.

The speculum will collect a lot of fluid so use that for inoculation as well

Inoculate and incubate

Always working inside the candle star, take the cooled petri dishes (check that they are actually hardened).

Using an inoculator, pick up the sample you are interested in (from a swab, speculum, underwear, etc...) and smear it on the surface of the medium drawing zigzag lines

Close the petri dish and seal with parafilm

Place in incubator upside down, setting the incubator at 37°. Candida grows overnight while lactobacilli are slow and will take about three days to incubate. The growth cycles of the bacteria are about 12 hours so you can occasionally go and check how the colonies are growing.

These are the reference pictures we took to give you an idea of the difference between candida and lactobacilli. In some cases you can see both.

Candida

Lactobacilli

Maintenance of bacteria and medium

Bacteria grow when they have the necessary food and the ideal temperature. To stop their growth, you can put sealed petri dishes in the freezer by hibernating them. We recommend placing them in a polystyrene box and keeping them separate from food in any case.

Bacteria are not eternal, so if kept out of the freezer they will feed and form further colonies until they die. Unused medium will harden and you can keep it in the fridge for a few days. If you want to reuse it, put it in a microwave or water bath to fluidise it.

It is important to take care with their disposal. We suggest autoclaving the medium in order to deactivate bacterial growth and sterilise the materials. Place the petri dish (if glass) with the medium for 30 minutes at 121 degrees. Otherwise you can use alcohol solution. They are organisms so throw them in the wet and separate the medium from the plastic petri dish. If you use glass petri dishes, sterilise them in an autoclave and reuse them without any problems.

Note on responsibility and safety

Although extracted from your body, bacteria still retain their capacity for interaction and symbiosis with the ecosystem in which they are embedded and thus also with you and your space. We do not believe in a pathological view of our relationship with the micro-organisms that always surround us. The word Contamination is defined by the Italian dictionary as 'physical or moral contact that disrupts the hygienic balance or traditional or individual values; corruption'. The moral meaning of a biological interaction underlies a vision of a pure and positive superiority that must not be soiled. We do not believe in contamination as a moral injury but as a mixing of heterogeneous elements that can create something new that does not have a moral meaning but rather one of responsibility. Some bacteria and fungi that you grow in the petri dish or that contaminate (in the biological sense) the bacterial culture may interact with you and your body, or that of other people. We think it is important to reflect on this and take care of it. And why not open to further reflection on our responsibilities and care when our microbiomes interact with each other.

Here are some non-moral suggestions on how to take on this responsibility, but we would welcome further advice and thoughts:

- Take care of the space in which you work, cleaning objects and surfaces.
- Prepare special spaces and materials that you will use for these experiments and magic. And possibly separate them from what you use in your days in the kitchen or the rest of the space.
- The same applies to the use of the fridge and freezer. Don't mix bacteria and food.
- It is just a matter of balance and care; not scientific sterility.
- Bacteria are also subject to the law, which can be very restrictive on these experiments. For example, do not travel with bacteria :)
- Take care of those who may be in an immunological deficiency condition, i.e. whose immune system due to other diseases is exposed to high risks in case of contagion with viruses and bacteria.